

# **medTester 5000C** Automated Biomedical Equipment Test System

## **Technical Data**



The medTester 5000C is an automated system designed for electrical safety testing and performance verification. It is compatible with most Fluke Biomedical testing devices and a majority of the popular Computerized Maintenance Management Systems (CMMS) in the US. The med-Tester 5000C provides a completely integrated solution for standardized and streamlined testing and record keeping.

The medTester 5000C can accept inventory and equipment testing procedures from the CMMS and control Fluke Biomedical testing devices in the field. The tool automatically runs all the tests and stores detailed results. Technicians can either print data or transfer it to CMMS for simple, automatic record keeping and exceptional documentation. The medTester 5000C automated system can provide up to a 50-percent saving in time and money over manual testing.

Modules are available for each Fluke Biomedical testing device that the medTester 5000C controls, for expanded capability, according to specific needs.

The medTester 5000C also serves as a standalone device to run integrated automated or manual electrical safety tests. The tool comes with ten preprogrammed and five user-programmable electrical-safety-testing sequences for one-button testing of virtually every piece of equipment in the hospital.

### **Key features**

- Easy verification of biomedical equipment to manufacturer's specifications
- Ten preprogrammed and five user-programmable electrical-safety-testing sequences
- · Convenient transfer of equipment inventory and testing procedures from CMMS
- Module options to automate testing of most Fluke Biomedical testing devices
- · Automatic storage of detailed test results for printing or transfer to CMMS
- Compliant with ANSI/AAMI (1993) and NFPA-99 (2005) standards
- 20 A device testing with GFCI protection
- Wedge hardware option for extended serial port use, optional PC-style keyboard and barcode scan gun



## **Specifications**

Modes of operation				
Fully equipped, with four operational modes	Manual, autosequence, medCh	Manual, autosequence, medCheck, and remote control		
Input power supply	Line-voltage/frequency input	115 V ac ± 10%/60 Hz		
	Test-receptacle type	USA, 20 A		
System/line voltage				
Range (full scale)	200 V	200 V		
Accuracy	± 5% of range ± 1 LSD	$\pm$ 5 % of range $\pm$ 1 LSD		
Resolution	0.1 V	0.1 V		
Equipment current	· ·			
Range (full scale)	0 A to 20 A	0 A to 20 A		
Accuracy	± 5% of range	$\pm$ 5% of range		
Resolution	0.01 A	0.01 A		
Ground resistance	÷			
Range (full scale)	0 Ω to 2 Ω	0 Ω to 2 Ω		
Accuracy	± 1% of range	± 1% of range		
Resolution	0.001 Ω (1 mΩ)			
Current source	100 mA dc			
Measurement type	True four-terminal technique	True four-terminal technique		
Test leads	Kelvin (2) insulated clip	Kelvin (2) insulated clip		
Leakage-current/voltage grad	lient			
Ranges (full scale)	200 µA and 2000 µA or mV			
Accuracy	$\pm$ 2.5% of reading + 3 LSD from	$\pm$ 1 % of reading dc + 3 LSD from 48 Hz to 1 kHz $\pm$ 2.5 % of reading + 3 LSD from 1 KHz to 100 KHz $\pm$ 5 % + 3 LSD from 100 KHz to 1Hz		
Resolution	0.1 μA or 0.1 mV	0.1 μA or 0.1 mV		
Measurement type	True-rms (autoranging) (AC + DC or dc-only response)			
Test-load selection	ANSI/AAMI ES1 1993	ANSI/AAMI ES1 1993		
Test-load impedence	1000 $\Omega \pm 1$ % at dc	$1000 \ \Omega \pm 1 \%$ at dc		
Isolation test				
Test selection (full scale)	Patient leads to ground	Patient leads to ground		
Lead combinations	All leads; or individual leads – tested as a single lead)	All leads; or individual leads – RL, RA, LA, LL, and V1/V6 (V1 through V6		
Available current	Limited by internal 120 k $\Omega$ res	Limited by internal 120 k $\Omega$ resistor		
Resolution	0.1 μΑ	0.1 μΑ		
Ranges (full scale)	200 µA and 2000 µA	200 μA and 2000 μA		
ECG binding posts				

pins and disposable snap electrodes



Performance waveforms				
ECG performance test waves (lead I, Vp-p)	Square wave	2 Hz, 1 mV		
	DC pulse	4 s, 1 mV		
	Sine wave	0.5 Hz, 10 Hz, 40 Hz, 60 Hz, and 100 Hz, 1 mV		
	Square wave	1 kHz, 1 mV		
	Triangle	2 Hz, 1 mV		
	CMRR	60 Hz sine wave with 1 k $\Omega$ imbalance in LA		
	Normal sinus	30 BPM, 60 BPM, 120 BPM, and 240 BPM		
Arrhythmias				

Atrial fibrillation; second-degree A-V Block, Type 1; premature atrial contractions; missed beat at 80 BPM and 120 BPM; PVC 1 left; PVC 2 right; multifocal PVCs; PVC 1, R on T; A pair of PVCs; run of 5 PVCs; run of 11 PVCs, MF; right bundle branch block; ventricular tachycardia; ventricular fibrillation; asystole

Environmental requirements				
Operating temperature	15 °C to 55 °C (59 °F to 95 °F)			
Storage temperature	0 °C to 50 °C (32 °F to 122 °F)			
General information				
Clock/date functions	Time and date formats	24 hour (hh:mm:ss) and mm/dd/yy		
Display characteristics	Туре	80-character, alphanumeric liquid crystal display (LCD)		
	Size	2 lines x 40 characters		
	Backlight	LED with adjustable brightness control		
Dimensions (LxWxH)	25.4 cm x 35 cm x 10.2 cm (10 in x 13.8 in x 4 in)			
Weight	5 kg (11 lb)			



#### **Ordering information**

Models

5000C medTester 5000C (20 A, 115 V ac) 5000C/B medTester 5000C/B, CMMS Connectivity Bundle

#### medTester 5000C

**Standard accessories** 9508-0272 Operators Manual 9530-0045FG Soft Vinyl Accessory Pouch 9501-0032FG Two Kelvin Cables 9503-0004FG Two Ground-Pin Adapters

#### medTester 5000C/B

**Standard accessories** 

**5000C** medTester 5000C 5000C-02 RS-232/Printer Module 5000C-03 100-Record Storage Module 5000C-04 Expanded Record Storage Module 5000C-06 Data-Transfer Module 5000C-07 medCheck Module 5000C-10 CMMS Interface Module

#### **Optional accessories**

5000C-02 Performance Enhancement Module 2: RS-232 Printer 5000C-03 Performance Enhancement Module 3: 100 Records

5000C-04 Performance Enhancement Module 4: Expanded Memory 5000C-05 Performance Enhancement Module 5: Waves/ **Extended Test** 5000C-06 Performance Enhancement Module 6: Data Transfer 5000C-07 Performance Enhancement Module 7: medCheck 5000C-08 Performance Enhancement Module 8: Defibrillator 5000C-09 Performance Enhancement Module 9: IV Pump 5000C-10 Performance Enhancement Module 10: CMMS Interface 5000C-11 Performance Enhancement Module 11: ESU 5000C-12 Performance Enhancement Module 12: Sp02 5000C-13 Performance Enhancement Module 13: Pacer 5000C-14 Performance Enhancement Module 14: NIBP 5000C-WEDGE Wedge Adapter (eight 25 in serial ports, as well as AT or PS/2 keyboard port) 9513-0212 Mini PC-Style External Keyboard (83 keys, AT or PS/2, wedge adapter required 9513-0221 Laser Barcode Gun (wedge adapter required) 5000C-PRINTER 5000C-PRINTER, Brady TLS Test Label Printer Kit medTester 5000C V 5.10 or greater and 115 V ac only 9530-0066 Multipurpose Hard-sided Carrying Case for medTester 5000C without wedge adapter Call Interface Cables for specific test-device connection

#### About Fluke Biomedical

Fluke Biomedical is the world's leading manufacturer of quality biomedical test and simulation products. In addition, Fluke Biomedical provides the latest medical imaging and oncology quality-assurance solutions for regulatory compliance. Highly credentialed and equipped with a NVLAP Lab Code 200566-0 accredited laboratory, Fluke Biomedical also offers the best

in quality and customer service for all your equipment calibration needs. Today, biomedical personnel must meet the increasing regulatory pressures, higher quality standards, and rapid technological growth, Biomedical provides a diverse range of software and hardware tools to meet today's challenges.

Fluke Biomedical Regulatory Commitment As a medical test device manufacturer, we recognize and follow certain quality standards and certifications when developing our products. We are ISO 9001 and ISO 13485 medical device certified and our products are:
NIST Traceable and Calibrated
UL, CSA, ETL Certified, where required

#### Fluke Biomedical.

Better products. More choices. One company.

**FLUKE** 

Biomedical

#### Fluke Biomedical 6045 Cochran Road Cleveland, OH 44139-3303 U.S.A.

Fluke Biomedical Europe Science Park Eindhoven 5110 5692EC Son, The Netherlands

For more information, contact us: In the U.S.A. (800) 850-4608 or Fax (440) 349-2307 In Europe/M-East/Africa +31 40 267 5435 or Fax +31 40 267 5436 From other countries +1 (440) 248-9300 or Fax +1 (440) 349-2307 Email: sales@flukebiomedical.com Web access: www.flukebiomedical.com

©2006-2013 Fluke Biomedical. Specifications subject to change without notice. Printed in U.S.A. 12/2013 2804434E\_EN

Modification of this document is not permitted without written permission from Fluke Corporation.